

DISPOSABLE VIRUS SAMPLING TUBE


Disposable virus sampling tubes can collect , store and transport specimens to special people.



Model	Liquid Tubes	Reagents Volume	Non-inactivated	Inactivated
VTM-A-CM0101	2mL	1mL	<input checked="" type="checkbox"/>	
VTM-A-CM0202	5mL	2mL	<input checked="" type="checkbox"/>	
VTM-A-CM0203	5mL	3mL	<input checked="" type="checkbox"/>	
VTM-A-CM0403	10mL	3mL	<input checked="" type="checkbox"/>	
VTM-A-CM0405	10mL	5mL	<input checked="" type="checkbox"/>	
VTM-B-CM0101	2mL	1mL		<input checked="" type="checkbox"/>
VTM-B-CM0202	5mL	2mL		<input checked="" type="checkbox"/>
VTM-B-CM0203	5mL	3mL		<input checked="" type="checkbox"/>
VTM-B-CM0403	10mL	3mL		<input checked="" type="checkbox"/>
VTM-B-CM0405	10mL	5mL		<input checked="" type="checkbox"/>

Disposable virus sampling tube is an important in virus detection.

Our company produces inactivated and non inactivated virus sampling tubes for selection.



1ML

2ML

3ML

5ML

6ML

10ML

REAGENTS VOLUME

Application

