

Cell Medical (Dongguan) co, .Ltd

Viral Transport Medium Kits

Disposable virus sampling tubes can collect , store and transport specimens to special people.



Viral Transport Medium Kits

The VTM kits is available in 2 formats : Without and with swabs.

Specification

- A) 1pcs flocked swab ,1 pc 10ml tube with 3ml medium refilled
- B) 2pcs flocked swab ,1 pc 10ml tube with 3ml medium refilled
- C) 1pcs flocked swab ,1 pc 5ml tube with 2ml medium refilled
- D) 2pcs flocked swab ,1 pc 5ml tube with 2ml medium refilled

Nylon fiber tip, ABS plastic shaft



Throat flocked swab CM-93050,
152mm length,80mm break point



Tube Volume:15ml
Material PP



Tube Volume:10ml
Material PP



Nasal flocked swab
CM-96000,150mm length,
80mm break point



Tube Volume:5ml
Material PP



Tube Volume:2ml
Material PP




Mouth Flocked swab CM-98000,
147mm Length,30mm Break point

Model	Specification
VTM-B-CM040320S	2With Nasal Swab& Throat Swab ; 3ml inactivated medium ;10ml tube
VTM-B-CM040311S	With Nasal Swab ; 3ml inactivated medium ;10ml tube
VTM-B-CM040312S	With Throat Swab ; 3ml inactivated medium ;10ml tube
VTM-A-CM040320S	With Nasal Swab & Throat Swab ; 3ml Non inactivated medium ;10ml tube
VTM-A-CM040311S	With Nasal Swab ; 3ml Non- inactivated medium ;10ml tube
VTM-A-CM040311S	With Throat Swab ; 3ml Non- inactivated medium ;10ml tube

Note: 1ml, 2ml, 3ml, 5ml, 6ml, 10ml and other viral transport medium can meet different needs of users

Disposable virus sampling tube is an important in virus detection.

Our company produces inactivated and non inactivated virus sampling tubes for selection.



1ML

2ML

3ML

5ML

6ML

10ML

REAGENTS VOLUME

Application

