



# flocked Swabs

Application



# flocked Swabs

Medical nylon flocked specimen collection swab has the ability of superior specimen elution, it releases up to 90% of the specimen collected compared to 25% in a standard fiber tipped swab, this ensure high reliability of the results.

Cell Medical (Dongguan) co, .Ltd

## Oral Flocked Swab



CM-98000



Oral:CM-98000

D 4.8mm±0.2mm

Breaking point



# Oral Flocked Swab



CM-93050



## Swab size

Throat: CM-93050



## Nasal Flocked Swab




CM-96000



### Swab size

Nasal: CM-96000 80mm breakpoint



Nasal: CM-96000N 90mm breakpoint

